
PEM REVERSIBLE FUEL CELL C-7121

TECHNICAL SPECIFICATIONS OF THE CELL PEM

Function electrolyzer


Input voltage.....1,7V – 3V (DC)
Consumption.....0,7 A @ 2 V
Half of hydrogen production.....5ml/min.
average production of oxygen.....2,5 ml/min.

Fuel cell Function

Output voltage.....0,6 V (DC)
Output current.....300 mA.
Output power.....180 mW

CONTENT

- A PEM reversible cell
- B Black wire plugs
- C Red cable with plug
- D Two brass pins
- E Flexible tube
- F Two red caps
- G Two black color caps
- H Syringe for hydration

A

B

C

D

F

G

E

H

Size..... 54 x 54 x 17 mm

Weight.....14 gr.

INTRODUCTION

Cell PEM (Polymer Electrolyte Membrane) that combines the functions of electrolyzer and fuel cell into a single component.

Operation

reversible: when connected to a generator of electricity (photovoltaic solar cell, turbine generator with crank, batteries, etc ...) and oxygen to produce hydrogen from deionized water (acts as an electrolyzer). By applying an electric charge generated electricity from hydrogen (fuel cell acts like).

Perfect for science labs, technology classrooms and demonstrations on hydrogen technology. Also for power portable computers and robots.

GENERAL SAFETY INSTRUCTIONS

- 1 - This cell is intended only for adults who have enough knowledge of fuel cells to develop their own projects by themselves or over 14 supervised by an adult who has read and understood the instructions in this manual.
- 2 - Before using this equipment should read and understand these instructions.
- 3 - Tools are required to mount this unit. We must pay due attention when being handled tools to avoid injury.
- 4 - Some small pieces are fragile, need to be handled carefully to avoid breaking them.
- 5 - Do not use any portion, element or component of the cell supplied for any purpose other than those specified in this manual. The cell may not be disassembled.
- 6 - At the end of each use to drain all the components containing water, hydrogen and oxygen.
- 7 - Wash your hands after working or experimenting with fuel cells

WARNING FOR THE CORRECT USE OF THE CELL

- 1 - Use only distilled water. Using distilled water will not damage the electrodes of the fuel cell. The fuel cell uses a carbon-catalyst.
Platinum nano, these particles are very sensitive to impurities that are always present in the water not distilled. When you see that the fuel cell starts to rust, it means that their experiments did not use distilled water.
- 2 - Before making sure you've electrolysis cell well hydrated by injecting water into the syringe. Leave the water on the inside with 5-10 minutes for complete hydration.
- 3 - The first time you perform electrolysis may lagua falls short optimal results in the production of hydrogen and oxygen in ratio 2-1. The maximum yield was obtained after in the water repeated three or four times the entire electrolysis process. This is due to increased hydration of the PEM membrane of the cell after continued use.
Ideal working temperature: 20 ° C ~ 30 ° C. Before proceeding with electrolysis make sure deposits are filled with enough distilled water.
- 4 - The fuel cell must be hydrated using the search ONLY W side and NEVER through oped H. Failure to do so will result in blocking the flow of hydrogen
- 5 - Check nozzles d containers are not sealed by gas plastic chips on top. If the small nozzles are clogged, the accumulate too much pressure within the fuel cell and that harm it.
- 6 - When the fuel cell used many times, the water from the top of the tank might not fall to internal lcontenedor. This is due to cate has caused a vacuum in the tubes. Disconnect the tube upper nozzle the fuel cell and water fall to inside container
- 7 - The fuel cell outdoors is very sensitive to organic compounds volatile, which will affect their performance. When you being use the kit to avoid air intake is advisable to keep the cell within a resealable plastic bag (Ziploc ® type, Minigrip ®, etc., as the original packing bag). This will protect the fuel cell when not using it.
- 8 - When using a solar panel or other power source, the output current should not exceed 0.7 A and the normal stress should not be higher than 2V. A solar panel whose output voltage or current exceed destroy the fuel cell.
- 9 - Whatch when connecting cables to the fuel cell. If you reverse the wires completely destroy


SOLUTIONS A POSSIBLE PROBLEMS

- 1 - Water levels do not drop by disconnecting the tubes out of gas both sides of the fuel cell
SOLUTION : Check pipes and nozzles are not clogged
- 2- The electrolyzer n oproduce hydrogen and / or oxygen
SOLUTION : a) Check for loose connections and if the cables are correctly connected. The fuel cell can be damaged irreparably if the red cable the supply is connected to the black jack fuel cell.
b) Make sure the power switch circuit is connected (ON).
- 3 - The water electrolysis process is very slow
SOLUTION : a) Using the syringe, inject distilled water into the oxygen side of the fuel cell. Wait 5 minutes before starting the process of electrolysis.
b) If using batteries, replace them with new ones.
c) If you use a solar Panel, check if it is properly oriented and enough sun shines.
d) If you use a wind turbine make sure that enough wind blows
- 4 - The application stops working, but there is hydrogen in the tank
SOLUTION : a) Purge gas and perform the electrolysis of water for 4-5 minutes.
Disconnect the tubes to purge oxygen and hydrogen gases. Make a new electrolysis until the container is full. Now connect the application to Cell combustible.If the problem persists, go to the next step.
b) Leave the electrolysis process last about 10 minutes, until reduced the remaining water. Purge gas to expel and the water of fuel cell.
Make a new electrolysis of water until the container is full. Then connect the application to fuel cell.

NOTE: This kit is recommended for children aged 14 years, always accompanied by an adult

Cebekit[®] is a registered trademark of the group Fadisel

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